

## Webinar Series in TEM: Introduction to Transmission Electron Microscopy Part 1

#### Riza Iskandar

Application Scientist

Material Science

Material and Structural Analysis - APAC



Direktorat Riset & PENGEMBANGAN



The world leader in serving science

## **Various Microscopes**

Microscopy is one of the few methodologies applied to nearly every field of science and technology in use today

**Thermo Fisher** 

SCIENT



# **Introduction to Transmission Electron Microscopy**

Outline: Part 01

- TEM Basics
- SEM vs TEM: What are the differences?
- Conventional TEM Imaging
  - Bright-Field Imaging
  - Dark-Field Imaging
- Electron Diffraction
  - Selective Area Electron Diffraction
  - Convergence Beam Electron Diffraction
  - Nano Beam Electron Diffraction
  - High-resolution TEM

## **Introduction to Transmission Electron Microscopy**

#### Outline: Part 02

- Scanning Transmission Electron Microscopy (STEM)
- Spectroscopy in TEM
- Tomography in TEM: For 2D to 3D Imaging
- Sample Image Analysis
  - Good and Reliable TEM Data: How to produce and Analyze
  - Example of TEM Data
- TEM Facility in Universitas Indonesia

## **TEST QUIZ**

## Thermo Fisher

## Which of the following is a figure of TEM?







С

thermoscientific

А



## Which one is the image taken from TEM?







Α

В

С

# **Introduction to Transmission Electron Microscopy**

Outline: Part 01

- TEM Basics
- SEM vs TEM: What is differences
- Conventional TEM Imaging
  - Bright-Field Imaging
  - Dark-Field Imaging
- Electron Diffraction
  - Selective Area Electron Diffraction
  - Convergence Beam Electron Diffraction
  - Nano Beam Electron Diffraction
- High-resolution TEM

Electron Microscope History

- 1897: Thompson describes the existence of negatively charged particles (electrons)
- 1925: De Broglie theorized that electrons have wave-like characteristics, addressing the wave/particle duality
- 1927: Thompson and Reid demonstrated the wave nature of electrons by diffraction experiments
- 1931: Ruska et al. build the first electron microscope (Nobel Prize in 1986)





### What is TEM?



٠

- To see **small objects** which cannot be seen with naked eyes, light microscope or even a SEM.
- To obtain structural information of **small objects**
- To analyze the chemical compositions of small objects



**Electron – Sample Interactions** 

#### **Thermo Fisher** s c | e N T | F | C

## Thermo Fisher

## **TEM Basics**

#### What can TEM Do?





Electron Diffraction pattern (SAED) taken from nano particles (left) as shown in TEM-BF (right)



Cat bone marrow with Feline Immunodeficiency Virus



EDX mapping of Cu and Fe elements distribution in steel





HRSTEM and HREDS SrTiO<sub>3</sub>





Au-Ni-Cl nanoparticles

Magnetic Hexa-ferrite

#### How does TEM Work?

#### Slide Projector vs TEM



Slide Projector	Transmission Electron Microscope
Visible light	electron
Glass lenses	Electron-magnetic lenses
$\lambda = 450-650 \text{ nm}$	$\lambda$ = 0.0025 nm (200 kV)
spatial resolution = 100 nm $(d > \lambda)$	spatial resolution = 0.24 nm $(d \cong \lambda)$

## Why Use Electrons?

- The resolution of light microscopy is limited by the wavelength of visible light (400 700nm)\*
- Electrons, that are particles as well as a wave, have much shorter wavelength, which gives much better resolution
- De Broglie equation:  $\lambda = h/mv$

Wave	U (kV)	Relativistic (λ=nm)	$r_{\rm th} = 1.22\lambda/\beta$
	100 120	0.0037	δ
$\overline{}$	200	0.0025	
$ \rightarrow \lambda \rightarrow$	300	0.0020	image plane

\*X-ray wavelength is about 0.05-0.25 nm

## **Electron Sources**

Thermionic



#### **Field Emission Gun (FEG)**





 $\text{Tungsten}(W)/\text{LaB}_6 \rightarrow \text{SFEG} \rightarrow \text{XFEG} \rightarrow \text{X-CFEG}$ 

Lifetime SFEG is 5 times longer than  $LaB_6$ LaB<sub>6</sub> = 1000 hours vs SFEG = 5000 hours Brighter sources Narrower energy spread Smaller emitter tip Better Performance



## **Optics: Image and Diffraction Formation**



- Object plane: a plane contains object points
- Back focal plane (BFP): the focal plane lies behind the lens
- Image plane: a plane contains image points





## QUIZ

- 1. A TEM is used to:
  - A. See small objects that cannot be seen with other imaging tools
  - B. Obtain structural information of small objects
  - c. Analyze the chemical compositions of small objects
  - D. A, B and C
- 2. A good TEM performance will depend on:
  - A. A higher acceleration voltage
  - B. Better electron sources
  - c. An expert operator
  - D. **A, B, C**
  - E. A and B

## QUIZ

- 3. What from these statements are NOT true for FEG electron sources:
  - A. It has a smaller size of emitter tip than Thermionic electron sources
  - B. It has a longer operating time than Thermionic electron sources
  - c. It has a higher brightness level than Thermionic electron sources
  - D. It has a broader energy spread than Thermionic electron sources
- 4. The most important lens in the TEM is:
  - A. Condenser lens
  - B. Objective lens
  - c. Intermediate lens
  - D. Projective lens

# **Introduction to Transmission Electron Microscopy**

Outline: Part 01

- TEM Basics
- SEM vs TEM: What are the differences?
- Conventional TEM Imaging
  - Bright-Field Imaging
  - Dark-Field Imaging
- Electron Diffraction
  - Selective Area Electron Diffraction
  - Convergence Beam Electron Diffraction
  - Nano Beam Electron Diffraction
- High-resolution TEM

#### What are the differences?

#### SEM: Scanning Electron Microscope



## **TEM:** Transmission Electron Microscope





#### What are the differences?

## SEM: Scanning Electron Microscope

- Smaller/shorter
- Acceleration voltage: 5kV 30 kV
- Resolution  $\geq$  0.7 nm
- Focused scanning beam
- Larger specimen chamber
- Larger samples



## **Transmission Electron Microscope**

- Larger/taller
- Acceleration voltage: 60-300kV or 30kV 1MV\*
- Resolution ≤ 0.1 nm
- Broad static beam and focused scanning beam (STEM)\*
- Smaller specimen chamber
- Thin samples of  $\approx 100 \text{ nm}$



## Samples







- Diameter: 2-3 mm
- Thickness: 100 150 nm

## QUIZ

- 1. Which statements are true?
  - A. TEM's sample is more difficult to prepare than SEM's sample
  - B. Most SEM detectors are located below the sample
  - c. One can achieve resolution below 1 nm with SEM
  - D. Scanning mode is only used on SEM
- 2. Which of these are significant differences between SEM and S/TEM imaging?
  - A. Sample preparation
  - B. Image resolution
  - c. Contrast mechanism
  - D. Accelerate voltage

# **Introduction to Transmission Electron Microscopy**

Outline: Part 01

- TEM Basics
- SEM vs TEM: What are the differences?
- Conventional TEM Imaging
  - Bright-Field Imaging
  - Dark-Field Imaging
- Electron Diffraction
  - Selective Area Electron Diffraction
  - Convergence Beam Electron Diffraction
  - Nano Beam Electron Diffraction
- High-resolution TEM

## **TEM Bright Field Image**

SAED Patterns and respective TEM-BF with and without Objective Aperture (OA)







Edge of Objective Aperture



**Thermo Fisher** 

SCIENTI

# **TEM Bright Field Image**

Bright Field Imaging Formation

- Bright field is the standard mode for imaging
- Selected by the objective aperture (TEM) or the detector (STEM)
- In TEM mode, on-axis dark-field (beam tilting) is better than offaxis dark-field (aperture shifting)

Objective aperture (OA): an aperture that is used to produce TEM-BF or TEM-DF







# **Bright and Dark Field Imaging**

**Bright Field Experiment Important Factors** 

- Sample is in Eucentric height
- Use the thinner part of the sample
- Refocus the image for fine-tuning (Gaussian Focus)
- Switch to diffraction mode
- If the sample is single crystal, try to orient the sample into a zone axis
- Insert and center the objective aperture
- A smaller objective aperture will give more contrast
- The image can be acquired

- Eucentric height: a sample position where the center of the sample image does not shift during specimen tilt
- Gaussian focus: a focus (imaging) condition for an ideal lens without aberrations
- Zone axis: a direction to a group planes that parallel to certain direction

## **TEM-Dark Field Image**

## Bright and Dark Field Imaging

TEM-BF

100 nm



**ThermoFisher** SCIENTIFIC

## **TEM-Dark Field Image**

Dark Field Imaging Formation

- Standard mode for imaging without the transmitted beam
- Loss of resolution due to higher C<sub>s</sub> at off-axis positions
- Two types of DF-Images: Off-axis and On-axis DF







Thermo Fisher



# Why TEM-DF is Important?

#### Benefit of Dark Field Image

- Grain Orientation
- Defect Analysis
- Phase formation





(a) Planes near an edge dislocation bend into the orientation for diffraction (b) BF image and (c) DF image of dislocations under a two-beam condition in an AI thin film. The inset in (b) shows the SAED pattern indicating the orientation condition for BF imaging.

## QUIZ

1. Which one of the following images is the real TEM-BF?



**Thermo Fisher** S C I E N T I F I C

- A. |
- в. II
- C. III
- D. **I,II**
- E. II, III

## QUIZ





## QUIZ

- **3.** With TEM-DF image, one can not obtain information of:
  - A. Crystal Defects
  - B. Grain Orientation
  - c. Element Identification
  - D. Phase Distribution
- **4.** Which is valid for acquiring a TEM BF:
  - A. The sample must be in eucentric height
  - B. The sample must be in zone axis
  - c. The smallest objective aperture must be used
  - D. The camera shortest camera length must be used to increase the contrast



# **Introduction to Transmission Electron Microscopy**

Outline: Part 01

- TEM Basics
- SEM vs TEM: What are the differences?
- Conventional TEM Imaging
  - Bright-Field Imaging
  - Dark-Field Imaging
- Electron Diffraction
  - Selective Area Electron Diffraction
  - Convergence Beam Electron Diffraction
  - Nano Beam Electron Diffraction
- High-resolution TEM

## **Electron Diffraction**

## X-rays vs Electron Diffraction

X-rays vs Electron Diffraction				
Parameter	X-rays	Electron		
Analysis	Intensity and position	Position		
Type of wavelength for single crystal	Multi-wavelength	Single-wavelength		
Acquisition time	Minutes to hours	Less than second		



Important facts about electron for diffraction:

- Electrons have shorter  $\lambda$  than X-rays
- Electron scattered more strongly
- Electron beams are easily directed







Single crystalline SAED pattern

SAED pattern Poly-crystalline SAED pattern

Amorphous SAED pattern

# **Type of Electron Diffraction in TEM**

#### Three Main Electron Diffraction Techniques



Selected Area Electron Diffraction (SAED)



Convergence Beam Electron Diffraction (CBED)





#### SAED Procedure



- Make sure sample in eucentric height
- As the name, the aperture is always needed to select an area and produce the Electron Diffraction Pattern
- Always take above three images when acquiring SAED pattern

Orienting a single crystal using the diffraction pattern: Single Crystal



• A double tilt holder is needed!

Far from zone axis





**Final adjusted SAED Pattern** 

**Diffraction and Interference** 



I hermo F

General concept: monochromatic & coherent radiation

- elastic scattering
- constructive interference only in certain directions
- special case:  $q_1 = 90^\circ$ : high resolution TEM

Bragg Equation

## 2 $d \sin \theta = n \lambda$

hermo

*d* is the spacing of atoms which scatter the electrons*n* is the order of diffraction*θ* the Bragg angle

For electron diffraction, electron wavelength is small,  $\theta$  is therefore small.

We expect diffraction from planes of atoms almost parallel to the incident electron beam.

(An Overview) SAED Pattern Analysis



For small angle of diffraction:  $r/L = 2\theta$ But  $\lambda = 2 d\theta$  $r/L = \lambda/d$  $r d = L \lambda$ 

Thermo Fisher

#### $L \lambda$ is the camera constant

#### SAED Pattern Analysis Procedures



- Make sure that the microscope is calibrated
- Measure the distance from reflections or rings to the main spot (d/2)

Thermo Fi

- Inversed the measured value
- Compare the data (dspacing) with database to find the phase presence

# **Convergence beam electron diffraction (CBED)**

#### More Than Just Electron Diffraction

SAED pattern of Si <111>



- Small area  $\leq 100 \text{ nm}$
- CBED gives quantitative data on
  - Specimen thickness
  - Crystallographic data such as unit cells, Bravais lattice, crystal system and 3D full symmetry
  - Precise lattice-strain measurements
  - Valence-electron distribution, structure factor, and chemical bonding
  - Characterization online and planar defects

CBED pattern of Si <111>



Ideal for relatively thicker sample

# **Nano Beam Electron Diffraction**

#### Electron Diffraction on Small Volume



- Nanocrystals of Fe<sub>3</sub>O<sub>4</sub> (magnetite), which have been incorporated in melt spinning polyvinylidene fluoride (PVDF) fibers.
- Experimental NBED in STEM mode (a) with corresponding simulated pattern (b).
- The frame colors refer to the position of electron beam during acquisition of the diffraction patterns.
- Domains or particles can be analyzed at nm-range by collection electron diffraction pattern.

Data countersy of N. Wirch Central Facility for Electron Microscopy, RWTH Aachen University

## **Electron Diffraction**

## QUIZ

- 1. Which is not the properties of electron for diffraction:
  - A. Electron has shorter  $\lambda$  than X-rays
  - **B.** Electron is less destructive than X-rays
  - c. Electron scattered more strongly
  - D. Electron beams are easily directed
- 2. The conditions that are needed to analyze SAED patterns:
  - A. The SAED pattern has been calibrated
  - B. The aperture must be inserted
  - c. The crystal must be oriented in the zone axis
  - D. The sample must be in eucentric height

## **Electron Diffraction**

## QUIZ

- 3. Convergence beam electron diffraction (CBED) can be used to produce diffraction on the size of:
  - A. 500 nm
  - в. 100 nm
  - c. 1 nm
- 4. What is true about CBED patterns?
  - A. Use an aperture to select an area of interest
  - **B.** Can give information about specimen thickness
  - c. Can provide information about chemical bounding
  - D. It is more suitable for a thicker TEM sample

# **Introduction to Transmission Electron Microscopy**

Outline: Part 01

- TEM Basics
- SEM vs TEM: What are the differences?
- Conventional TEM Imaging
  - Bright-Field Imaging
  - Dark-Field Imaging
- Electron Diffraction
  - Selective Area Electron Diffraction
  - Convergence Beam Electron Diffraction
  - Nano Beam Electron Diffraction
- High-resolution TEM

Electron as a Wave

- Electron beam is a wave with *amplitude* (A) and phase (φ)
- The periodicity is the wave-length (e.g. 0.0025nm at 200kV) or, in terms of phase,  $2\pi$



**Imaging Formation** 



Parallel Incoming electron beam (wave)

High-resolution TEM: Contrast Transfer Function (I)

- Point resolution is the function of Cs and  $\lambda$
- It is the same for W, LaB6 and FEG
- But the information limit is different



- Point Resolution or Resolution:

   A wave number at which the phaseconstant transfer function at
   Scherzer focus intersects at phasezero axis
- Information limit:

Indicates the wave number at which phase information carried by the phase-contrast transfer function disappears.

High-resolution TEM: Contrast Transfer Function (II)

• CTF greatly depends on the focus value (Cs=1.2mm)



**Thermo Fisher** 

# **Advanced TEM Imaging: High-resolution TEM**

High-resolution TEM: What do we See on a HRTEM image?

Usually, you can not say where the atom is, but you can tell the distance on atomic scale and crystal defects

Thermo Físhei



## HRTEM and HR-STEM Imaging



HRTEM and HR-STEM images showing the atomic structure of Si<sub>3</sub>Ni<sub>4</sub>

• Si+4 N-3

## **Thermo Fisher**

# **High-resolution TEM**

## A good and Reliable HRTEM Image

- Sample:
  - Very thin  $\leq 50 \text{ nm}$
  - Free from contamination
- Sample position:
  - Eucentric height
  - Oriented in one low zone axis ([001], [120])
- Instrument:
  - Stable enough, i.e. HT on not less than 2 hours
  - No or less sample drift
  - The beam is well aligned
  - Highest magnification as possible\*
  - Adequate acquisition time
  - Focus series applied for acquisition



## QUIZ

- 1. To produce High-resolution TEM images (samples perspective):
  - A. Very thin
  - B. Crystalline
  - c. Oriented on a zone axis
  - D. Cannot work with powder samples
- 2. Point Resolution depends on the following factors:
  - A. Spherical Chromatic Aberration (Cs)
  - B. Electron Sources
  - c. Wavelength ( $\lambda$ )
  - D. Operation Voltage



## QUIZ

- **3.** The information that we can get from HRTEM images:
  - A. The exact atomic position
  - B. The structure of the sample
  - c. The type of atom
  - D. The distance on an atomic scale and crystal defects

# Thank you

56 mohamadriza.iskandar@thermofisher.com | 30-June-2021